

Clinical Measurements

# A Quick Guide to Capnography

Recording and analysis  
of the CO<sub>2</sub> waveform  
and its use in  
differential diagnosis

Prof. B. Smalhout  
M.D., Ph.D.

**PHILIPS**

---

*This guide provides convenient views of a selection of CO<sub>2</sub> waveforms, along with interpretation to explain their relationship to other routine physiological monitoring waveforms.*

*These illustrated waveforms are semi-schematic diagrams presented in their ideal shape.*

---

# Contents

## The CO<sub>2</sub> waveform

Analyzing the CO <sub>2</sub> waveform	2
Important basic rules	4

## The Capnogram

The normal capnogram	12
Expiratory problems	13
Curare capnogram	14
Cardiogenic oscillations	15
Camel capnogram	16
Iceberg capnogram	17
Leakage in the respiratory system	18
High end-tidal CO <sub>2</sub>	19
Low end-tidal CO <sub>2</sub>	23
Some disturbances of the respiratory rhythm	26
Patient fights the ventilator	32

## Differential diagnosis

Seven examples with interpretations	35
-------------------------------------	----

### Analyzing the CO<sub>2</sub> waveform

The CO<sub>2</sub> waveform can be analyzed for five characteristics:

- Height
- Frequency
- Rhythm
- Baseline
- Shape


The normal end-tidal value is approximately: 38 mmHg or 5% etPCO<sub>2</sub>

- Height depends on the end-tidal CO<sub>2</sub> value etPCO<sub>2</sub>
- Frequency depends on the respiratory rate
- Rhythm depends on the state of the respiratory center or on the function of the ventilator
- Baseline should be zero
- There is only one normal shape (see page 12)

## The CO<sub>2</sub> waveform

The waveforms are typically recorded or displayed at two different speeds:

***Real-time (high) speed at 12.5 mm/sec***



***Trend (slow) speed at 25 mm/min***




## Important basic rules

- 1. A sudden drop** in  $\text{CO}_2$  to zero or to a low level always indicates a technical disturbance or defect:
  - Spontaneous breathing or ventilated patients
    - Kinked ET-tube
    - $\text{CO}_2$  analyzer defective
  - Ventilated patients
    - Total disconnection
    - Ventilator defective
  
- 2. A sudden change** in baseline, sometimes combined with changes in plateau level, indicates:
  - Calibration error
  - $\text{CO}_2$  absorber saturated (check capnograph with room air)
  - Water drops in analyzer or condensation in airway adapter


**CHECK CAPNOGRAPH !**

# Important basic rules

1.



2.



## Important basic rules

- 3. A sudden decrease in  $\text{CO}_2$  value** (but not to zero) with spontaneous or ventilated breathing indicates:
  - Leakage in the respiratory system (low airway pressure)
  - Obstruction (high airway pressure)
- 4. An exponential decrease in  $\text{CO}_2$**  (washout curve) within one or two minutes always indicates a sudden disturbance in lung circulation or ventilation:
  - Circulatory arrest
  - Embolism
  - Sudden decrease in blood pressure
  - Sudden severe hyperventilation




## Important basic rules

3.



4.




## Important basic rules

- 5. Gradual increase in  $\text{CO}_2$**  for spontaneous or controlled breathing indicates:
  - Developing hypoventilation
  - Absorption of  $\text{CO}_2$  from peritoneal cavity (laparoscopy)
  - Rapidly rising body temperature
- 6. Sudden increase in  $\text{CO}_2$**  (spontaneous or controlled ventilation):
  - Injection of sodium bicarbonate
  - Sudden release of tourniquet (legs, arms, etc.)
  - Sudden increase in blood pressure (e.g., intravenous adrenaline)

## Important basic rules

5.



6.




## Important basic rules

- 7. Gradual upshift** in the  $\text{CO}_2$  baseline and topline can result from:
- Saturation of  $\text{CO}_2$  absorber
  - Calibration error
  - Technical error in  $\text{CO}_2$  analyzer
  - Increasing dead-space, resulting in re-breathing
  - $\text{CO}_2$  absorber switched off
- 8. Gradual lowering** of the end-tidal  $\text{CO}_2$ . The curve retains its normal shape but the height of the plateau falls gradually. In an artificially ventilated patient, this can be caused by:
- Gradual hyperventilation
  - Lowering body temperature
  - Decreasing body or lung perfusion

## Important basic rules


7.



8.



### The normal capnogram




There is only one normal shape.

Characteristics:

- Rapid increase from P to Q
- Nearly horizontal plateau between Q and R (slightly sloping up to R)
- From R rapid decrease to zero
- Points P, Q, R, and S appear as rounded corners. (P, Q, R is the expiratory phase. R, S, P is the inspiratory phase.)
- Slope of the plateau depends on the condition of the airways and lung tissue
- End-tidal value is **only** equivalent to the alveolar CO<sub>2</sub> when a nearly horizontal plateau is seen

## Representative capnograms

### Expiratory problems



In general, any airway obstruction limiting expiration.

Possible explanations:

- Kinked tube (developing out of a previously normal shape)
- Foreign body
- Herniated cuff (developing out of a previously normal curve)
- Bronchospasm
- Emphysema
- Bronchial asthma

## Representative capnograms

### “Curare” capnogram\*




- Caused by lack of coordination between intercostal muscles and diaphragm
- Note the cleft in right third of the plateau. CO<sub>2</sub> mostly too high. The depth of the cleft is proportional to the degree of muscle paralysis. Seen only in spontaneous respiration, or when the patient starts to fight the ventilator
- Also seen in patients with cervical transverse lesions

\* The expression “curare capnogram” was given at the time that curare was a generally accepted muscle relaxant (about 1960). Nowadays, other muscle relaxants are in use, but the shape of the capnographic curve in patients who are partially paralyzed is still the same.



## Representative capnograms

### Cardiogenic oscillations




Caused by the beating of the heart against the lungs:

- Small tidal volume in combination with low respiratory rate
- At the end of a very long expiration

Can be caused by a central depression of the respiratory system or by the ventilator running too slowly.

## Representative capnograms

### Camel capnogram



- Can be seen in patients in the lateral position on the operating table
- During either spontaneous or controlled respiration

## Representative capnograms

### Iceberg capnogram




Caused by a combination of a muscle relaxant and a central acting analgesic drug (e.g., morphine, fentanyl, etc.):

- Mixture of cardiogenic oscillations and “curare” cleft (see page 14)
- No plateau. Low respiratory rate
- CO<sub>2</sub> higher than normal
- Seen only in spontaneous respiration

## Representative capnograms

### Leakage in the respiratory system




Irregularity mostly in expiratory limb (see arrows):

- Shape and site of this irregularity depend on the localization and severity of the leak in the anesthetic system (cuff, valves, tubing, etc.). CO<sub>2</sub> could be too high due to hypoventilation or too low due to the addition of leaking air
- Other possible leakage shapes have to be differentiated from other disturbances (e.g., camel curve (see page 16), etc.)

## High end-tidal CO<sub>2</sub>

### A. With normal respiratory rate




Normal plateau but higher than normal end-tidal CO<sub>2</sub>. Can be seen in artificially ventilated patients:

- When the ventilator is running at normal respiratory rate but the minute volume is too low
- With primary normal respiratory rate and minute volume but with a rapidly rising body temperature (e.g., in malignant hyperthermia)

## High end-tidal CO<sub>2</sub>

### B. With bradypnea




Long plateau but higher than normal CO<sub>2</sub>. Respiratory depression without an attempt to compensate can be seen:

- In cases of high ICP or respiratory depression due to morphinic-based drugs (e.g., pethidine, fentanyl, etc.)
- In artificially ventilated patients, when the ventilator is running with both respiratory rate and minute volume too low

## High end-tidal CO<sub>2</sub>


### C. Hypoventilation with tachypnea



Short plateau but higher than normal end-tidal CO<sub>2</sub>. Respiratory depression with an attempt to compensate by higher respiratory rate can be caused by:

- Volatile anesthetics during spontaneous respiration (e.g., halothane, etc.)
- Ventilator running at high rate but with low tidal volume

### D. Very severe hypoventilation




Very severe hypoventilation (very low tidal volume) with a high respiratory rate as an attempt to compensate. A misleading low level of CO<sub>2</sub> is recorded on the capnogram. Mostly no proper plateau. After thorax compression or forced exhalation (see arrow), true CO<sub>2</sub> value becomes visible:

- Can be seen under spontaneous or controlled breathing
- In patients with spontaneous respiration but with severe respiratory paralysis caused by paralyzed respiratory muscles
- Malfunction of ventilator or leakage in respiratory system



## Low end-tidal CO<sub>2</sub>

### A. With normal respiratory rate




With normal respiratory rate and plateau but a lower than normal end-tidal CO<sub>2</sub>. Can be observed in artificially ventilated patients:

- When the ventilator runs with a normal rate but the minute volume is too high
- Who are in shock
- With normal respiratory rate and tidal volume but with a low body temperature

Can also be seen in patients with spontaneous respiration when they are compensating a metabolic acidosis.

## Low end-tidal CO<sub>2</sub>

### B. With bradypnea




With bradypnea but a lower than normal end-tidal CO<sub>2</sub> and a long plateau in:

- Artificially ventilated patient when the ventilator is running at a low rate and a high minute volume
- In patients with spontaneous respiration with:
  - Damage to the central nervous system (e.g., the so-called central neurogenic hyperventilation)
  - A low body temperature and respiratory depression caused by analgesics

## Low end-tidal CO<sub>2</sub>

### C. With tachypnea



With tachypnea but a lower than normal end-tidal CO<sub>2</sub> and short plateau:

- In patients on artificial ventilation when the ventilator is running at high rate and with a high minute volume
- In patients with spontaneous respiration who are:
  - In pain
  - Trying to compensate a metabolic acidosis
  - Hypoxic
- In some cases of central neurogenic hyperventilation
- In severe shock conditions

## Some disturbances of the respiratory rhythm

### Cheyne-Stokes respiration



Only seen with spontaneous respiration. Cardiac oscillations (indicated by the arrows) after every respiratory group.

- Can be seen in cases of severe cerebral arteriosclerosis, brain damage, intoxication, etc.

## Some disturbances of the respiratory rhythm

### Heaving respiration




Can be seen in some patients as a transition condition between Cheyne-Stokes and normal breathing.

- Only seen with spontaneous respiration. The tidal volume changes regularly, hence the waving characteristic of the upper limit on the capnographic trend recording.

## Some disturbances of the respiratory rhythm

### Gasping respiration




- Very low respiration rate (2-6/min)
- CO<sub>2</sub> mostly higher than normal
- Often cardiogenic oscillations after every capnographic curve (see arrows)

Seen in very severe respiratory depression or in dying patients

## Some disturbances of the respiratory rhythm

### **Very irregular or chaotic respiration**




- No regularity
- All curves differ in size, shape, and height. Average CO<sub>2</sub> level above normal

Seen in severe cerebral damage

## Some disturbances of the respiratory rhythm

### Sighing respiration





- A regular pattern with regular intervals interrupted by a deep sigh (indicated by arrows)
- Physiological in babies, small children, and very old people during sleep or when under an anesthetic
- Pathological in young people and adults when deep sighs are more frequent than once in 5 minutes. An indication of brain damage
- Average  $\text{CO}_2$  level can be normal, high, or lower than normal



## Some disturbances of the respiratory rhythm

Can also be seen in ventilated patients when the ventilator has an intermittent deep sigh mechanism. In normal lungs, the deep sigh  $\text{CO}_2$  level is lower than average (A). In cases of obstructive lung disease, the deep sigh level will be higher than average (B).




## Patient fights the ventilator

When a patient starts to breathe against the ventilator (A), the regular pattern of the capnogram is interrupted. The respiratory activity of the patient increases quickly. The end-tidal  $\text{CO}_2$  rises slightly due to the increasing metabolism of the contracting respiratory muscles.

Capnogram created by the ventilator (B).

Capnogram created through the attempted spontaneous respiration of the patient (C).

- During an anesthetic this capnogram indicates that another dose of muscle relaxant should be given to the patient



---

**Note:**

Combinations of the waveforms discussed up to this point are always possible. Results depend on the clinical condition of the patient and the technical status of the instrumentation used.

The iceberg capnogram is an example of a combination waveform. (see page 17)



### **Differential diagnosis**

Correct interpretation of the many possible capnographic curves can only be achieved by comparison with other parameters recorded simultaneously. This makes a differential diagnosis possible.

Useful parameters for this purpose are:

- ECG/heart rate
- Blood pressure (direct or indirect)
- Body temperature
- Plethysmogram (taken from the earlobe or finger)
- PaCO<sub>2</sub>
- PaO<sub>2</sub>
- Airway pressure
- Central venous pressure
- Acid-base status


## **Example 1**

- ECG: disappears, no heart rate
- Pleth: becomes a straight line
- Blood pressure: drops
- Capnogram: no change

### ***Interpretation:***

Technical defect when the vital signs monitor (ECG, pleth, blood pressure) and the capnograph are separate instruments.

Patient is not in danger. Without proper circulation a normal capnogram is not possible.



## Example 2

- Controlled respiration, temperature 37°C
- ECG: normal
- Pleth: normal
- Blood pressure: constant
- Capnogram: rapid decrease in CO<sub>2</sub>

### **Possible interpretations:**

- Considerable leakage in the respiration system. **PATIENT IN DANGER**
- Technical disturbance in capnograph. **Patient not in danger**

### **Differential diagnosis:**

- Check airway pressure
- Check capnograph with own expiratory air (always around 38 mmHg or 5%)



ECG



ECG: normal

PLETHYSMOGRAM



Plethysmogram: normal

BLOOD PRESS



Blood pressure: unchanged

CAPNOGRAM



Capnogram: shows rapid decrease

12.5 mm/sec 25 mm/min

### Example 3

- ECG: increased heart rate, occasional PVCs
- Pleth: diminishing amplitude and irregularities in amplitude
- Blood pressure: increasing
- Capnogram: CO<sub>2</sub> level rises

### **Possible interpretations:**

- Adrenaline intoxication (e.g., local anesthesia)
- Effect from manipulating a pheochromocytoma
- Effect from painful stimulus
- Patient awakes
- **PATIENT POSSIBLY IN DANGER**

ECG



ECG: increased heart rate, occasional PVCs

PLETHYSMOGRAM



Plethysmogram: diminished amplitude and irregularities in amplitude

BLOOD PRESS



Blood pressure: increasing

CAPNOGRAM



Capnogram: CO<sub>2</sub> rises

12.5 mm/sec 25 mm/min

### **Example 4**

- ECG: asystole, after some PVCs
- Pleth: diminishing amplitude developing into a straight line
- Blood pressure: drops to zero
- Capnogram: washout curve. Value drops towards zero but holds at a level of several mmHg

### **Interpretation:**

This combination is typical for cardiac arrest even though there can still be some electrical activity in the ECG.

**PATIENT IN MORTAL DANGER**

ECG



ECG: asystolia often some PVCs

PLETHYSMOGRAM



Plethysmogram: amplitude/height amplitude and becomes straight line

BLOOD PRESS



Blood pressure: drops to zero

CAPNOGRAM



Capnogram: wash-out curve. Drops mostly not completely to zero.  
Remains some time on a level of several mm's Hg

12.5 mm/sec 25 mm/min

## Example 5

- ECG: tachycardia, ST depression. Or bradycardia with very low blood pressure
- Pleth: decrease in amplitude
- Blood pressure: decreases
- Capnogram: drops in parallel with blood pressure

### **Possible interpretations:**

- Severe blood loss
- Severe circulatory collapse through other causes:
  - Anaphylactic shock
  - Cardiac malfunction
  - Overdose of certain cardio-depressing drugs (e.g., halothane or barbiturates)
- **PATIENT IN DANGER**

ECG



ECG: tachycardia, ST depression. Or bradycardia with very low blood pressure

PLETHYSMOGRAM



Plethysmogram: decrease of amplitude

BLOOD PRESS



Blood pressure: decreases

CAPNOGRAM



Capnogram: drops parallel with blood pressure

12.5 mm/sec 25 mm/min

## Example 6

- ECG: no changes. Sometimes bradycardia develops
- Pleth: rapid increase in amplitude
- Blood pressure: rapid decrease, remains low
- Capnogram: no significant change

### **Possible interpretations:**

- Effect of neuroleptic drugs (e.g., Droperidol), alpha-blocker or ganglion blockers or anesthetic drugs
- Induction of anesthesia

Blood pressure decreases due to vasodilatation. Perfusion improved.

Patient is not in danger as long as the capnogram remains unchanged.



EKG



EKG: no changes.  
Sometimes bradycardia develops

PLETHYSMOGRAM



Plethysmogram: rapid increase  
in amplitude

BLOOD PRESS



Blood pressure: rapid decrease,  
remains low

CAPNOGRAM



Capnogram: no significant change.

12.5 mm/sec

25 mm/min

## Example 7

- ECG: remains unchanged at first although hypoxic changes may occur after a few minutes with PVCs
- Pleth: develops into an almost straight line (often within ten seconds)
- This may be preceded by a broadening of the record for a few seconds
- Blood pressure: falls within thirty seconds (often to a very low level)
- Capnogram: falls within one minute (often to a very low level of  $\text{etCO}_2$ )

### **Interpretation:**

Large pulmonary embolism (often air embolism). Recovery is usually gradual. Even after a small air embolism, which is not life-threatening, the capnogram does not return to its original level for at least five to ten minutes.

**PATIENT IN SEVERE DANGER**

ECG



ECG: regular, often  
some PVCs

PLETHYSMOGRAM



Plethysmogram: diminished  
amplitude and becomes  
straight line

BLOOD PRESS



Blood pressure: drops to  
low level

CAPNOGRAM



Capnogram: wash-out curve.  
Drops mostly not completely  
to zero. Remains on a level  
of several mm's Hg

12.5 mm/sec 25 mm/min

## **Author:**

Prof. B. Smalhout, M.D., Ph.D.  
Professor of Anesthesiology  
University Hospital Utrecht  
The Netherlands

## **Other Publications**

- Atlas of Capnography  
Smalhout/Kalenda  
Kerckebosch - Zeist  
The Netherlands
- The Suffocating Child  
Smalhout/Vaughan  
Boehringer-Ingelheim  
Germany



**Philips Medical  
Systems is part  
of Royal Philips  
Electronics**

www.medical.philips.com  
medical@philips.com  
fax: +31 40 27 64 887

Philips Medical Systems  
3000 Minuteman Road  
Andover, MA 01810-  
1085  
(800) 934-7372

© Koninklijke Philips  
Electronics N.V. 2004  
All rights are reserved.  
Reproduction in whole  
or in part is prohibited  
without the prior  
written consent of the  
copyright holder.

Philips Medical Systems  
Nederland B.V. reserves  
the right to make  
changes in specifications  
and/or to discontinue any  
product at any time  
without notice or  
obligation and will not be  
liable for any  
consequences resulting  
from the use of this  
publication.

4522 981 93901/862  
SEP 2004